

iOSBluetoothSDK V1.3_Release

Date: 0723, 2020

Author: Liang Fang

SDK Version: 1.3

MCU: 2.3

Update record:

1. V1.3 fix the bug of failed callback in iOS13.5 system
2. V1.2 supports multiple or single connections
3. V1.0 supports single connection

Content

iOSBluetoothSDK V1.3_Release	1
iOSBluetoothSDK Development Guide	4
Introduction.....	4
Your First Project: iOSBluetoothSDK V1.3_Release Demo.....	5
iOSBluetoothSDK V1.3 API Reference.....	10
HZLBlueData Reference	10
Blue3OrBlue4ManagerReference	12

iOSBluetoothSDK Development Guide

Introduction

This guide will teach you how to use HZLBlueTooth SDK for iOS to write iOS applications that can acquire brainwave data from MacroTelligence 's Hardware (BrainLink Pro & BrainLink Lite) . This will enable your iOS apps to receive and use brainwave data such as BLEMIND and BLEGRAVITY acquired via Bluetooth, MacroTelligence 's Hardware and File source encapsulated as HZLBlueTooth. HZLBlueTooth SDK for iOS supports upgrading Hardware

Function:

Receive brainwave data. One or more Bluetooth devices can be connected at the same time.

Files included:

- API Reference (this document)
- SDK static library and headers
- iOSBluetoothSDK V1.3_Release.a
- HZLBlueData.h
- Blue3OrBlue4Manager.h
- iOSBluetoothSDK V1.3_Release Demo

Supported devices:

- Data format with power
 - BrainLink_Pro
 - Jii
- Data format without power
 - BrainLink_Lite
 - Mind Link


iOS VerrSION :

- iOS 9.0 +

Your First Project: iOSBluetoothSDK V1.3_Release Demo

Step 1:


1.1 Import the iOS framework libraries CoreBluetooth.framework and ExternalAccessory.framework in the Build Phases of TARGETS in the Xcode project:



Add com.neurosky.thinkgear In the Info.plist: (ios13 needs to add the Bluetooth permission privacy - Bluetooth always usage description and privacy - Bluetooth peripheral usage description)

Key	Type	Value
▼ Information Property List		
Localization native development region	String	\$(DEVELOPMENT_LANGUA
Executable file	String	\$(EXECUTABLE_NAME)
Bundle identifier	String	\$(PRODUCT_BUNDLE_IDEN
InfoDictionary version	String	6.0
Bundle name	String	\$(PRODUCT_NAME)
Bundle OS Type code	String	APPL
Bundle versions string, short	String	1.0
Bundle version	String	1
Application requires iPhone environment	Boolean	YES
▶ Required background modes	Array	(3 items)
Launch screen interface file base name	String	LaunchScreen
Main storyboard file base name	String	Main
▶ Required device capabilities	Array	(1 item)
▼ Supported external accessory protocols	Array	(1 item)
Item 0	String	com.neurosky.thinkgear
▶ Supported interface orientations	Array	(3 items)
▶ Supported interface orientations (iPad)	Array	(4 items)

1.2 If you want Bluetooth to work in the background, please set it as follows ,
Don't set it if you don't need it:



Step 2:

Import header file

```
#import "HZLBlueData.h"
```

```
#import "Blue3OrBlue4Manager.h"
```

Function one : Receive data

```
// Bluetooth multi-connection
```

```
// blueNames: connectable device
```

```
NSArray *blueNames = @[@"BrainLink",@"BrainLink_Pro",@"jii@jii-***"];
```

```
[Blue3OrBlue4Manager logEnable:YES];
```

```
[[Blue3OrBlue4Manager sharedInstance] configureBlueNames:blueNames ableDeviceSum:7];
```

```
// Connect bluetooth callback successfully
```

```
__weak FactoryViewController *weakSelf = self;
```

```
[Blue3OrBlue4Manager sharedInstance].blueConBlock = ^(NSString * markKey){
```

```
    // Determine connected devices 123456 corresponds to ABCDEF respectively
```

```
    if ([markKey isEqualToString:@"1"]){
```

```
        NSLog(@"Bluetooth device A is connected");
```

```
    }
```

```
    .....
```

```
    else if ([markKey isEqualToString:@"6"]){
```

```
        NSLog(@"Bluetooth device F is connected");
```

```
    }
```

```
};
```

```
// Bluetooth disconnect callback
```

```
[Blue3OrBlue4Manager sharedInstance].blueDisBlock = ^(NSString * markKey){
```

```
    // Determine disConnected devices
```

```
    if ([markKey isEqualToString:@"1"]){
```

```
        NSLog(@"Bluetooth device A is disConnected");
```

```
    }
```

```

.....
    else if ([markKey isEqualToString:@"6"]){
        NSLog(@"Bluetooth device F is disConnected");
    }
}
};

```

//蓝牙数据回调

```

[Blue3OrBlue4Manager sharedInstance].hzlblueDataBlock_A = ^(HZLBlueData *blueData,
BlueType conBT, BOOL isFalseCon) {
    if (conBT == BlueType_Pro) {
        if (blueData.bleDataType == BLEMIND) {
            weakSelf.ALabl.text = [NSString stringWithFormat:@"sigal:%d att:%d med:%d
ele:%d ap:%d del:%d theta:%d lowAlp:%d highAlp:%d lowBe:%d highBe:%d lowGa:%d highGa:%d version:%d
grid=%d",blueData.signal,blueData.attention,blueData.meditation,blueData.batteryCapacity,blueData.
ap,blueData.delta,blueData.theta,blueData.lowAlpha,blueData.highAlpha,blueData.lowBeta,blueData.
highBeta,blueData.lowGamma,blueData.highGamma,blueData.hardwareVersion,blueData.grind];
            // when the signal value is 0, the bluetooth device is worn
            // note: if the bluetooth device is connected but not worn, Greater than 0 and less
than or equal to 200
            if(blueData.signal == 0){
                weakSelf.ASignalIV.image = [UIImage
imageNamed:@"signal_zhengChang"];
            }else{
                weakSelf.ASignalIV.image = [UIImage imageNamed:@"signal3.png"];
            }
        }
        else if (blueData.bleDataType == BLEGRAVITY) {
            weakSelf.ACircleLabl.text = [NSString stringWithFormat:@"x:%d y:%d
z:%d",blueData.xvlaue,blueData.yvlaue,blueData.zvlaue];
        }
        else if (blueData.bleDataType == BLERaw) {
            weakSelf.ARawLabl.text = [NSString stringWithFormat:@"raw:%d
eye:%d",blueData.raw,blueData.blinkeye];
        }
    }
    else if (conBT == BlueType_Jii){
        if (blueData.bleDataType == BLEMIND) {
            weakSelf.ALabl.text = [NSString stringWithFormat:@"sigal:%d att:%d med:%d
ele:%d
ap:%d",blueData.signal,blueData.attention,blueData.meditation,blueData.batteryCapacity,blueData.ap
];
            if(blueData.signal == 0){
                weakSelf.ASignalIV.image = [UIImage
imageNamed:@"signal_zhengChang"];
            }else{

```

```

        weakSelf.ASignalIV.image = [UIImage imageNamed:@"signal3.png"];
    }
}
}
else if (conBT == BlueType_Lite) {
    if (blueData.bleDataType == BLEMIND) {
        weakSelf.ALabl.text = [NSString stringWithFormat:@"sigal:%d att:%d
med:%d del:%d theta:%d lowAlp:%d highAlp:%d lowBe:%d highBe:%d lowGa:%d
highGa:%d",blueData.signal,blueData.attention,blueData.meditation,blueData.delta,blueData.theta,blu
eData.lowAlpha,blueData.highAlpha,blueData.lowBeta,blueData.highBeta,blueData.lowGamma,blueD
ata.highGamma];

        //信号值为 0 即佩戴了蓝牙设备
        //注：如果连接了蓝牙设备而未佩戴，信号值为大于 0 且小于或等于 200
        if(blueData.signal == 0){
            weakSelf.ASignalIV.image = [UIImage
imageNamed:@"signal_zhengChang"];
        }else{
            weakSelf.ASignalIV.image = [UIImage imageNamed:@"signal3.png"];
        }
    }
    else if (blueData.bleDataType == BLERaw) {
        weakSelf.ARawLabl.text = [NSString stringWithFormat:@"raw:%d
eye:%d",blueData.raw,blueData.blinkeye];
    }
}

if (isFalseCon) {
    NSLog(@"A device has a false connection");
}

};
.....

[Blue3OrBlue4Manager sharedInstance].hzlblueDataBlock_G = ^(HZLBlueData *blueData, BlueType
conBT, BOOL isFalseCon) {
    .....
};

[[Blue3OrBlue4Manager sharedInstance] connectBlue3OrBlue4];

// Active bluetooth disconnect
[[Blue3OrBlue4Manager sharedInstance]disconnectBlue3OrBlue4];

```


iOSBluetoothSDK V1.3 API Reference

HZLBlueData Reference

Overview

The HZLBlueData class is a data model

Enum

```
typedef enum : NSUInteger {
    BlueType_NO = 0,
    BlueType_Lite,
    /*The current connection is the BrainLink_Lite data format device with BLEMIND and BLERaw type data*/
    BlueType_Pro,
    /* The current connection is the BrainLink_Pro data format device with BLEMIND, BLEGRAVITY, and BLERaw type data*/
    BlueType_Jii,
    /* The current connection is Jii*/
}BlueType;

typedef NS_ENUM(NSUInteger,BLEDATATYPE){
    BLEMIND    =    0,           // basic brain wave data
    BLEGRAVITY,           // gravity data
    BLERaw,           // blink data
};
```

Basic Brainwave Data :

- signal,
- attention,
- meditation,
- delta,
- theta,
- lowAlpha,
- highAlpha,
- lowBeta,
- highBeta,
- lowGamma,
- highGamma,
- ap,
- batteryCapacity,
- hardwareVersion,
- grind

Gravity Sensor Data:

- xvlaue,
- yvlaue,
- zvlaue

Raw& Blink Data:

- raw,
- blinkeye

Note :

When Jii is connected, only signal, attention, meditation, batteryCapacity, ap data type is available.

When BrainLink_Lite is connected, only signal , attention , meditation , delta , theta , lowAlpha , highAlpha , lowBeta , highBeta , lowGamma , highGamma , raw , blinkeye data type is available.

Instructions of some Instance Property

- **signal:** It represents the signal value of the MacroTelligence 's Hardware. When the signal is 0, it means that the MacroTelligence 's Hardware has been put on, and when the signal is greater than 0 and less than or equal to 200, it means that the MacroTelligence 's Hardware is connected to the iPhone.
- **batteryCapacity:** In percentage terms. minimum value is 0, maximum value is 100
- **ap:** Appreciation value
- **hardwareVersion:** Hardware version. The first version value is 255 , when you update the MacroTelligence 's Hardware , the version value will be smaller.
- **xvlaue :** gravity value in The x axis (Pitching Angle)
- **yvlaue :** gravity value in The y axis (Yaw Angle)
- **zvlaue :** gravity value in The z axis (Roll Angle)

Blue3OrBlue4ManagerReference

Overview

The Blue3OrBlue4Manager class handles interaction between a MacroTelligence's Hardware and an iOS device.

Instance Property

Successful callback of bluetooth connection

```
@property (nonatomic,copy)Blue3OrBlue4Connect blueConBlock;
```

Bluetooth disconnect callback

```
@property (nonatomic,copy) BlueConnectdismiss blueDisBlock;
```

Note: devices are connected in the order of A B C D E F. Six data callbacks (hzluedatablock_A ...) are used to ensure the independence of data. The data between various devices can be accepted at the same time without mutual influence. Bluetooth 4.0 devices can connect up to six, but it is difficult to connect successfully.

If you want to use a single connection, the input parameter of abledevicesum is 1. Only call hzluedatablock_A.

Data callback for every device.

```
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_A;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_B;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_C;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_D;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_E;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_F;
```

Connection status of every device

```
@property (nonatomic,assign)BOOL connected_A;  
@property (nonatomic,assign)BOOL connected_B;  
@property (nonatomic,assign)BOOL connected_C;  
@property (nonatomic,assign)BOOL connected_D;  
@property (nonatomic,assign)BOOL connected_E;  
@property (nonatomic,assign)BOOL connected_F;
```

Method

Print log does not print by default

```
+ (void)logEnable:(BOOL)enable;
```

Initialization (singleton)

```
+ (instancetype)shareInstance;
```

Parameter Configuration :

Parameter interpretation :

blueNames: **Device name (Bluetooth 4.0 device) or MFI (Bluetooth 3.0 device) that can be connected**

```
NSArray *blueNames = @[@"BrainLink",@"BrainLink_Pro",@"jii@jii-***"];
```

1.jii@jii-*** : indicates the device name with the prefix of jii - that can be connected , 'jii@' means jii device. The content After @ is the device name. '***' indicates the prefix is the same.

2.BrainLink : Indicates the device name or MFI called brainlink that can be connected.

ableDeviceSum: Number of Bluetooth devices that can be connected

```
-(void)configureBlueNames:(NSArray *)blueNames ableDeviceSum:(int)deviceSum;
```

Connect bluetooth device

```
-(void)connectBlue3OrBlue4;
```

Disconnect bluetooth device

```
-(void)disConnectBlue3OrBlue4;
```

Manual test of false connection (definition of false connection: when signal is equal to 0 and the continuous 10 values of attention and mediation remain unchanged, it is considered to be a false connection. SDK will disconnect the Bluetooth connection of the current device and automatically connect again)

```
-(void)testAFalseCon:(BOOL)isTest; //Manual test A device false connection
```

```
-(void)setTestToZero;//Cancel all manual test false connections
```