# HZLBluetooth\_V1.2 SDK for iOS

---

Date: 0421, 2020

Author: Liang Fang

SDK Version: 1.2

MCU: 2.3

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---

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
```
#import "HZLBlueData.h"
```

```
#import "Blue3OrBlue4Manager.h"
```

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```
NSArray *blueNames = @[@"BrainLink",@"BrainLink_Pro",@"jii@jii-***"];
[Blue3OrBlue4Manager sharedInstance].blueConBlock = YES;
[[Blue3OrBlue4Manager sharedInstance] configureBlueNames:blueNames ableDeviceSum:6];
```

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```
__weak FactoryViewController *weakSelf = self;
[Blue3OrBlue4Manager sharedInstance].blueConBlock = ^( NSString *markKey){
    //判断连接的设备 123456 分别对应 ABCDEF
    if ([markKey isEqualToString:@"1"]){
        NSLog(@"A 设备蓝牙连接成功");
    }
    .....
    else if ([markKey isEqualToString:@"6"]){
        NSLog(@"F 设备蓝牙连接成功");
    }
};
```

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```
[Blue3OrBlue4Manager sharedInstance].blueDisBlock = ^( NSString *markKey){
    //判断断开的设备
    if ([markKey isEqualToString:@"1"]){
        NSLog(@"A 设备蓝牙断开");
    }
};
```

```

    }
    .....
    else if ([markKey isEqualToString:@"6"]){
        NSLog(@"F 设备蓝牙断开");
    }
}
};

//蓝牙数据回调
[Blue3OrBlue4Manager sharedInstance].hzlblueDataBlock_A = ^(HZZBlueData *blueData,
BlueType conBT, BOOL isFalseCon) {
    if (conBT == BlueType_Pro) {
        if (blueData.bleDataType == BLEMIND) {
            weakSelf.ALabl.text = [NSString stringWithFormat:@"sigal:%d att:%d med:%d
ele:%d ap:%d del:%d theta:%d lowAlp:%d highAlp:%d lowBe:%d highBe:%d lowGa:%d highGa:%d version:%d
grid=%d",blueData.signal,blueData.attention,blueData.meditation,blueData.batteryCapacity,blueData.
ap,blueData.delta,blueData.theta,blueData.lowAlpha,blueData.highAlpha,blueData.lowBeta,blueData.
highBeta,blueData.lowGamma,blueData.highGamma,blueData.hardwareVersion,blueData.grind];

            //信号值为 0 即佩戴了蓝牙设备
            //注：如果连接了蓝牙设备而未佩戴，信号值为大于 0 且小于或等于 200
            if(blueData.signal == 0){
                weakSelf.ASignalIV.image = [UIImage
imageNamed:@"signal_zhengChang"];
            }else{
                weakSelf.ASignalIV.image = [UIImage imageNamed:@"signal3.png"];
            }
        }
        else if (blueData.bleDataType == BLEGRAVITY) {
            weakSelf.ACirlceLabl.text = [NSString stringWithFormat:@"x:%d y:%d
z:%d",blueData.xvlaue,blueData.yvlaue,blueData.zvlaue];
        }
        else if (blueData.bleDataType == BLERaw) {
            weakSelf.ARawLabl.text = [NSString stringWithFormat:@"raw:%d
eye:%d",blueData.raw,blueData.blinkeye];
        }
    }
    else if (conBT == BlueType_Jii){
        if (blueData.bleDataType == BLEMIND) {
            weakSelf.ALabl.text = [NSString stringWithFormat:@"sigal:%d att:%d med:%d
ele:%d
ap:%d",blueData.signal,blueData.attention,blueData.meditation,blueData.batteryCapacity,blueData.ap
];

            //信号值为 0 即佩戴了蓝牙设备
            //注：如果连接了蓝牙设备而未佩戴，信号值为大于 0 且小于或等于 200
            if(blueData.signal == 0){
                weakSelf.ASignalIV.image = [UIImage

```

```

imageNamed:@"signal_zhengChang"];
        }else{
            weakSelf.ASignalIV.image = [UIImage imageNamed:@"signal3.png"];
        }
    }
}
else if (conBT == BlueType_Lite) {
    if (blueData.bleDataType == BLEMIND) {
        weakSelf.ALabl.text = [NSString stringWithFormat:@"sigal:%d att:%d
med:%d del:%d theta:%d lowAlp:%d highAlp:%d lowBe:%d highBe:%d lowGa:%d
highGa:%d",blueData.signal,blueData.attention,blueData.meditation,blueData.delta,blueData.t
heta,blueData.lowAlpha,blueData.highAlpha,blueData.lowBeta,blueData.highBeta,blueData.lo
wGamma,blueData.highGamma];
        //信号值为 0 即佩戴了蓝牙设备
        //注：如果连接了蓝牙设备而未佩戴，信号值为大于 0 且小于或等于
200
        if(blueData.signal == 0){
            weakSelf.ASignalIV.image = [UIImage
imageNamed:@"signal_zhengChang"];
        }else{
            weakSelf.ASignalIV.image = [UIImage
imageNamed:@"signal3.png"];
        }
    }
    else if (blueData.bleDataType == BLERaw) {
        weakSelf.ARawLabl.text = [NSString stringWithFormat:@"raw:%d
eye:%d",blueData.raw,blueData.blinkeye];
    }
}

if (isFalseCon) {
    NSLog(@"A 设备假连接");
}
};
.....

[Blue3OrBlue4Manager sharedInstance].hzlblueDataBlock_F = ^(HZLBlueData *blueData, BlueType
conBT, BOOL isFalseCon) {
    .....
};

[[Blue3OrBlue4Manager sharedInstance] connectBlue3OrBlue4];

// 主动断开蓝牙
[[Blue3OrBlue4Manager sharedInstance]disconnectBlue3OrBlue4];

```

# HZLBluetooth API 参考

---

## HZLBlueData 参考

### Overview

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### Enum

```
typedef enum : NSUInteger {
    BlueType_NO = 0,
    BlueType_Lite,
    /*连接的是 BrainLink_Lite 数据格式设备 ,有 BLEMIND、BLERaw 类型数据 */
    BlueType_Pro,
    /*连接的是 BrainLink_Pro 数据格式设备 ,有 BLEMIND、BLEGRAVITY、BLERaw 类型数据 */
    BlueType_Jii,
    /*连接的是 Jii*/
}BlueType;

typedef NS_ENUM(NSUInteger,BLEDATATYPE){
    BLEMIND    =    0,           //脑波数据
    BLEGRAVITY,           //重力数据
    BLERaw,           //Raw 眨眼数据
};
```

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- signal, 设备佩戴质量
- attention, 专注度
- meditation, 放松度
- delta,
- theta,
- lowAlpha,
- highAlpha,
- lowBeta,
- highBeta,
- lowGamma,
- highGamma,
- ap, 喜好度
- batteryCapacity, 电池电量百分比

- hardwareVersion, 设备固件版本
- grind

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- xvlaue,
- yvlaue,
- zvlaue

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- raw,
- blinkeye

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连接 Jii , 只有 signal , attention , meditation , batteryCapacity , ap

连接 BrainLink\_Lite , 只有 signal , attention , meditation , delta , theta , lowAlpha , highAlpha , lowBeta , highBeta , lowGamma , highGamma , raw , blinkeye

**Instructions of some Instance Property**

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### Instance Property

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```
@property (nonatomic,copy)Blue3OrBlue4Connect blueConBlock;
```

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```
@property (nonatomic,copy) BlueConnectdismiss blueDisBlock;
```

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```
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_A;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_B;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_C;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_D;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_E;  
@property(nonatomic,copy)Blue3OrBlue4DataBlock hzlblueDataBlock_F;
```

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```
@property (nonatomic,assign)BOOL connected_A;  
@property (nonatomic,assign)BOOL connected_B;  
@property (nonatomic,assign)BOOL connected_C;  
@property (nonatomic,assign)BOOL connected_D;  
@property (nonatomic,assign)BOOL connected_E;  
@property (nonatomic,assign)BOOL connected_F;
```

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